

## Comparison Features of the MicroMed Electrosurgical Insulation Integrity Testers Models PD-8k and MM-5kV

<b>Instrument Model #</b>	<b>PD-8K</b>	<b>MM-5kV</b>
Portable/Battery Operated	Yes	Yes
Weight (pounds)	4.85	1.01
Ring Electrode for any size instrument	Yes	Yes
Brush Electrode	Yes	Yes
Reusable Ring /Brush	Yes	Yes
Test different electrosurgical instruments (e.g. laparoscopic and bi-polar)	Yes	Yes
On-going costs (*9 volt battery)	No	Yes*
Adjustable voltage to test various instruments jacket thicknesses	Yes	Yes
Multiple Alarms types (spark/ flash/ sound)	Yes	Yes
Adjustable Alarm sound	Yes	No
Easy to adjust knobs/Touch controls	Yes	Yes
UL Battery Adaptor	Yes	Yes (optional)
Battery	Special	9 volt
Battery pack (external)	No	Yes (optional)
Battery Life Indicator	Yes	Yes
Voltage Control	0-8kV	0-5kV
Dimensions (inches)	10.25 x 6.3 x 2.75	8.27 x 2.95 x 1.38
CE certified	Yes	Yes



[www.mcgan.com](http://www.mcgan.com)

**9 Liberty Hill Dr.  
Blackstone, MA 01504**

**Ph#: 508-876-1070 Fax #: 508-883-3434**